

# Magnetic Stirrer

## FEATURES


- COMPACT & LIGHTWEIGHT,
- Microprocessor Control
- Multi Direction ( Forward & Reverse )
- Knob / Soft Control ( Speed )
- SHORT CIRCUIT & OVERLOAD PROTECTION,
- Stainless Steel Top ( S.S )
- Dc ( External ) 12V - 2A , SMPS adapter
- Power Failure : Auto Restart set Speed



INJECTION MOLDING ABS PLASTIC CONSTRUCTION




BVS - MS1



Every minute it will run forward & reverse as per your set speed (Rpm)

BVS - MS2



( Memory at last set speed )

## Specifications

Model	BVS - MS1	BVS - MS2
Speed	1200 RPM	1200 RPM
Control	Knob type	Soft Touch Key
Capacity	2 Liter	2 Liter
Parameter	Speed & Direction	Speed & Direction
Motor	12 Volts - Dc	12 Volts - Dc
Soft Key	On / Off, Direction Control	Up , Down , On / Off, Direction
Multi Direction	Forward & Reverse	Forward & Reverse
Room Temperature	0 - 40°C	0 - 40°C
Related Humidity	0 - 90% non condensing	0 - 90% non condensing
Power (i/p)	Dc - 12V (2A )	Dc - 12V (2A )
Dimensions	156mm x 180mm x 40mm	156mm x 180mm x 40mm
Weight	900 gram	900 gram
Fuse	1A ⚡	1A ⚡

## OPTIONAL - ACCESSORIES

\* 12V - Dc motor

Dc power adapter ( SMPS )



Magnetic Bit's

