

Electrophoresis Power Supply


FEATURES

- COMPACT & LIGHTWEIGHT,
- 4 Digit Led Seven Segment,
- Knob Control & Membrane Key,
- 100 Hours Timer Control with Buzzer Indication,
- SHORT CIRCUIT & OVERLOAD PROTECTION,
- 36 MODELS DIGITAL POWER SUPPLY.



Injection Molding ABS Plastic Construction Panel Box

4 Digit Seven Segment - 100 Hours Selectable Timer
Real Time : Power Failure (Auto Restart Set Program)



This Power Supply of Modular design of Bio - electric focusing unit enables the user to carry out almost all of the flatbed electrophoresis techniques

| Specification | BVS - 100SS | BVS - 200SS | BVS - 500SS | BVS - 1SA |
|--------------------|----------------------------|----------------------------|----------------------------|----------------------------|
| VOLTAGE | 0 - 300 V ⚡ |
| CURRENT | 0 - 100 mA | 0 - 200 mA | 0 - 500 mA | 0 - 999 mA |
| TIMER | 1 min to 99hr 59min |
| DISPLAY | 4 Digit Seven Segment |
| RESOLUTION | 1 V, 1 mA, 1 min |
| OUTPUTS | 4 parallel (4 mm Socket) |
| PARAMETERS | Voltage, Current, Time | Voltage, Current, Time | Voltage, Current, Time | Voltage, Current, Time |
| Memory Option | NO (Option) | NO (Option) | NO (Option) | NO (Option) |
| PARAMETER CONTROL | Manual Mode | Manual Mode | Manual Mode | Manual Mode |
| COMPUTER CONTROL | NO (Option) | NO (Option) | NO (Option) | NO (Option) |
| WORK TEMPERATURE | 0 - 40°C | 0 - 40°C | 0 - 40°C | 0 - 40°C |
| RELATIVE HUMIDITY | 0 - 90%, non condensing |
| POWER REQUIREMENTS | 230 VAC, 50/60 Hz, |
| DIMENSIONS | 135mm x 190mm x 80mm |
| (W x L x H) WEIGHT | 1.2 kg | 1.2 kg | 1.2 kg | 1.2 kg |
| Fuse | 2 A (each) |

BVS ➤ 0 TO 60V - 3000mA Digital Power supply Also Available for Vertical gel unit

Mfd By

Biovision Services

P6 / E , 1st Cross, Opp BWSSB Water Tank,
3rd Stage, Peenya Industrial Area, Bangalore - 560058.
e-mail : biovision@gmx.com, info-biovision@gmx.com
Tel : 080 - 28364622, +91 - 95909 05557, 95909 06664
Website : [www.biovisionservices.com / in](http://www.biovisionservices.com/in)