


Electrophoresis Power Supply

FEATURES

- COMPACT & LIGHT WEIGHT
- LED Display , Digital Control
- Knob Control (Voltage)
- SHORT CIRCUIT & OVERLOAD PROTECTION,
- **36 MODELS DIGITAL POWER SUPPLY.**

Control - Method




Rocker Switch

Soft Key



5 Step Voltage Selectable - Rocker Switch Control



Model	BVS - Reg 1	BVS - Reg 2	BVS - Reg 3	BVS - Reg 5
Voltage	50, 100, 150, 200, 250 (5 Step)			
Current	100 mA	200 mA	300 mA	500 mA
Control	Rocker Switch			
Parameter	Voltage			
Dimensions	156mm x 180mm x 40mm			
Power	Single Phase - 230 V , 50 Hz			
Fuse	2A (Small Size)			
Weight	700 grms			

INJECTION MOLDING ABS PLASTIC CONSTRUCTION (Microprocessor Based)



Control - 5 Step Rocker Switch

This Power Supply of Modular design of Bio - electric focusing unit enables the user to carry out almost all of the flatbed electrophoresis techniques

6 Step Voltage Selectable - Soft Key Control



Model	BVS - Mini 1	BVS - Mini 2	BVS - Mini 3	BVS - Mini 5
Voltage	50, 100, 150, 200, 250, 300 V (6 Step)			
Current	100 mA	200 mA	300 mA	500 mA
Control	Soft Membrane Key			
Parameter	Voltage & Status (Start / Stop)			
Dimensions	156mm x 180mm x 40mm			
Power	Single Phase - 230 V , 50 Hz			
Fuse	2A (Small Size)			
Weight	700 grms			

INJECTION MOLDING ABS PLASTIC CONSTRUCTION (Microprocessor Based)



Control - 6 Step Soft Membrane Key

Mfd By

Biovision Services

P6 / E , 1st Cross, Opp BWSSB Water Tank,
3rd Stage, Peenya Industrial Area, Bangalore - 560058.
e-mail : biovision@gmx.com, info-biovision@gmx.com
Tel : 080 - 28364622, +91 - 95909 05557, 95909 06664
Website : [www.biovisionservices.com / in](http://www.biovisionservices.com/)



GO FOR INNOVATION