

Orbital Shaker


offers advance technology Orbital Lab Shakers are known for their durability and rotary swirling action used extensively in tissue culture work, Aeration of fermentations and various other chemical mixing procedures. BVS Series Orbital Shakers are built for continuous duty.

If you need a specific solution (Heavy Duty Large Volume Mixing) and OEM builds, [Biovision Services](#) has the ability to manufacture the solution.



BVS - OS1

- ▶ Microprocessor Control
- ▶ 40 to 200 / 250 / 300 / 400 Rpm
- ▶ Knob Control - Rpm Control
- ▶ No Timer Control
- ▶ Multi - Flask Holder
- ▶ Input AC - 230V, 50Hz



BVS - OS2

- ▶ Microprocessor Control
- ▶ 40 to 200 / 250 / 300 / 400 Rpm
- ▶ Knob Control - Rpm Control
- ▶ 4 Digit - Timer (99hr : 59min)
- ▶ Multi - Flask Holder
- ▶ Input AC - 230V, 50Hz

Ideal for growth of bacterial and yeast



BVS - OS3

◀ BVS - OS3

- ▶ **PID - RPM Controller**
(Microprocessor Control)
- ▶ Soft touch - Key
- ▶ 40 to 200 / 250 / 300 / 400 Rpm
- ▶ Continuous duty
- ▶ Multi - Flask Holder
- ▶ LCD - Jambo Display
- ▶ Timer (99hr : 59min)
- ▶ 20 Memory (R / W)
- ▶ Key : Rpm, Time, Up, Down, Start / Stop
- ▶ Input AC - 230V, 50Hz